



CardiLink

CardiLink® Features

Manage your VIVEST Power Beat Series AED remotely with CardiLink®.

Take control of your installed base. Make sure your AED is ready at all times.

Move towards a service-driven business model and create recurring revenues.





Access CardiLink® from your smartphone, your laptop or your PC.


Download the CardiLink app from the Apple AppStore or the Google PlayStore. Your service partner can invite you to the platform. Receive the invitation link via email, set and confirm your password and login to CardiLink®.



Track the realtime location of your VIVEST Power Beat Series AED 24/7 via GPS positioning.



The included global roaming SIM will connect to the strongest available cellular network.

Subscribe to the notifications in the CardiLink® platform. When the accelerometer detects movements, receive an SMS, e-mail or push notification. Press the link and follow the realtime location of your device on the map.



**Take control of your service.
Perform maintenance when it is
necessary.**

Rest assured and trust the system. CardiLink® will notify you when the self-test fails, pads expire or the battery is nearing the end.





Make sure to give others access to your lifesaving device. Locate the nearest AED. Saving a life is teamwork!

- 24/7 AED readiness on your phone
- Real location tracking via GPS
- Notifications for movements and events

Be ahead - with CardiLink®. The most comprehensive AED tracking and monitoring solution on the market.